

Technical Datasheet



Quick Release Thermic Lance Holder



SKU: 3LHOLDERQR

The brass body of the holder is fitted with a 3/8" BSP nipple for connection to the oxygen supply, a ball valve for control of oxygen, and a quick release carrier to which the required size probe for your lance is fitted.

| | |
|---------------------|------------------------------|
| Lance Compatibility | 8mm nominal bore (1/4" BSP) |
| | 10mm nominal bore (3/8" BSP) |
| | 20mm nominal bore (3/4" BSP) |

The appropriate lance is screwed into the threaded recess in the probe.

Issue No. 1

Distributing the tools to build stronger communities

T 02380 228668
 E sales@expressweldcare.co.uk
 W expressweldcare.co.uk

